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#### INTENDED USE

In vitro diagnostic medical device for performing the agglutination test in a test tube for the titration of antibodies according to Weil-Felix. Test results should always be interpreted in relation to the clinical context. For professional use only.

#### CLINICAL SIGNIFICANCE

Weil Felix test was developed at the beginning of '900 and is based on the revealing of antibodies against some species of Proteus that contain antigens cross-reacting with Rickettsia genus. Proteus vulgaris OX19 reacts with sera of people infected with Rickettsiae (typhus group e RMSF Rocky Mountains Spotted Fever), Proteus vulgaris OX2 reacts with sera of people infected with Rickettsiae (Spotted Fever Group) and Proteus mirabilis OXK reacts with sera of people with scrub typhus caused dalla Orientia tsutsugamushi.

#### PRINCIPLE OF THE METHOD

When a serum, containing the specific agglutinins, reacts with the homologous antigen, under optimized and controlled conditions, it produces a visible agglutination. This reaction can be performed in a test tube, on a slide or in a microplate. In the test, the serum is diluted by doubling in serology tubes. After the antigen is added, the reaction is incubated for the expected period of time and then read. The degree of agglutination depends on the concentration of antigen, antibody, salt composition and temperature.

#### COMPONENTS

All concentrations refer to ready to use reagent:

Suspension	Proteus OX19	1x20 ml
	Proteus OX2	1x20 ml
	Proteus OXK	1x20 ml

Unstained bacterial suspension at optimal concentration for tube test.

Contains Sodium azide 0,95 g/L.

\*Warning: Products that contain sodium azide may react with lead and copper in plumbing to form explosive metal azides. On disposal flush with large volumes of water to prevent azide build up.

#### Other material required, but not supplied

Serological tubes, Physiological saline, automatic pipettes.

#### Storage and Stability

 = Storage Temperature 2-8 °C

Store reagents and controls at 2-8° C avoiding direct light. Do not freeze. Under above conditions reagents are stable until expiration date indicated on the label. Slight variations in composition among batches will not affect test results.

#### PRECAUTIONS AND WARNINGS

1. Reagents and waste materials shall be disposed of in accordance with Community waste provisions or national or regional provisions.
2. Reagents may contain non-active components such as preservatives and detergents. The total concentration of these components is below the limits set out in Regulation 1272/2008 EC and subsequent amendments and additions.
3. It is recommended that the reagent be handled in accordance with the rules of good laboratory practice and that appropriate personal protective equipment be used.
4. Do not use the reagent if it is visibly degraded (e.g. presence of corpuscles).
5. All human samples shall be handled and disposed of as potentially infectious material.
6. The kit should only be used by qualified and properly trained technical personnel.
7. Diagnoses shall be carried out exclusively by authorised and qualified personnel.
8. Comply with national directives on occupational safety and quality assurance.
9. Use equipment that complies with current regulations.

#### Reporting of serious incidents

In the event of a serious incident in relation to the use of the device, please inform the manufacturer (via your distributor) and the competent authority of the European Union member state where the incident occurred. For other jurisdictions, reporting must be made in accordance with regulatory requirements. Reporting serious incidents helps provide more information about the safety of the diagnostic medical device.

#### REAGENT PREPARATION

Reagents are liquid ready to use. Avoid any chemical or bacteriological contamination. Mix well reagents before use. Suspension must appear uniform and free of visible particles.

#### SPECIMEN COLLECTION AND PRESERVATION

Not-inactivated sera must be used. Samples can be stored at 2-8° C for 48 h before testing. For longer time of storage samples must be frozen (-20° C). Highly hemolyzed, lipemic or contaminated samples cannot be used. Samples with the presence of fibrin should be centrifuged before testing.

#### PROCEDURE

##### Technique

- 1 -Bring reagents and materials to room temperature. Low temperature can reduce the method sensitivity.
- 2 -Swirl gently the reagent before using. Suspension must appear uniform and free from visible particles.
- 3- pre-dilute testing sera 1:25 with buffered PBS pH 7,2 and transfer 0.5 mL in tube 1.

Tube	MU	1	2	3	4	5	6	Control
PBS	mL	0.5	0.5	0.5	0.5	0.5	0.5	0.5
Serum(1:25)	mL	0.5						
Carry out two-fold dilutions transferring 0.5 mL from the tube n. 1 to the next, until tube N. 6, discarding last 0.5 mL								
Diagnostic suspension	mL	0.5	0.5	0.5	0.5	0.5	0.5	0.5
Mix and incubate a 37°C for 24 hours								
Serum dilutions		1:100	1:200	1:400	1:800	1:1600	1:3200	-

#### Results

**Positive Reaction (agglutination):** Sedimented cells with clear supernatant.

**Negative Reaction (absence of agglutination):** Uniformly cloudy suspension without sedimentation

#### Titer

The titer is defined as reciprocal dilution of serum showing a visible bacterial agglutination. Sometimes we can have a border-line effect; for example, a tube with less turbidity if compared with negative control showing , in the bottom, a sedimentation button. In this case the tube which shows the border-line effect, and not the previous tube, fix the titer.

#### Results interpretation

Generally, titers equal or higher than 1:100 are probative for the diagnosis of the disease.

If we consider + as a medium titer and ++ as a higher titer, the diagnosis can be made according to the following table:

Disease	Agent	Weil – Felix Reaction		
		OX 19	OX 2	OX K
Epidemic typhus	<i>R. prowazekii</i>	++	+	-
Murine typhus	<i>R. mooseri</i>	++	+	-
Rocky Mountains spotted Fever	<i>R. rickettsii</i>	+ / ++	++	- / +
Boutonneuse Fever	<i>R. conori</i>	+ / ++	++	-
Tsutsugamushi Fever	<i>R. orientalis</i>	-	-	++
Q fever	<i>R. burneti</i>	-	-	-

To have a correct diagnosis it is discriminating measuring a significant growth of the titer between two samples drawn in the 5th and 12th day from the beginning of the fever .

#### Validation Test

##### Specificity

Tests carried out with three different lots of bacterial suspension of Proteus OX19, OX2 and OX K, on samples of human negative serum have given repeatedly negative results.

##### Precision

Test of repeatability (within run) and reproducibility (between run) carried out with three different lots of bacterial suspension of Proteus OX19, OX2 and OX K on normal human samples and animal samples at known titre, have given constantly the expected results.

##### Sensibility

Tests carried out with three different lots of bacterial suspension OX19, OX2 ed OX K, on human positive sera, had confirmed the titer obtained with reference methods.



##### Quality Control

Agglutinable suspensions in the kit must be analyzed using a positive serum as control and PBS as Negative Control. If it is impossible to obtain positive and negative results respectively, suspensions or Control serum would be considered impaired.

##### Stability

Stability test, in real time, carried out on three different lots of each bacterial suspension, have confirmed the correct functionality of the reagent for at least 36 months to 2-8° C.

#### Symbols used in IFU and Packaging

 In vitro diagnostic medical device vitro	 Manufacturer
 Catalogue Number	 Instruction for use
 Lot Number	 Temperature limitation
 Expiration date	

#### Bibliografia

1. Castaneda M.R. - Bull WHO, 9, 399, 1953
2. Sonnerwirth A.C. - 1970, in Gradwohl's Clinical Laboratory Methods and Diagnosis 7th ed.,p. 1482, The CV Mosby Co., St. Louis

REVISION	DATE	CHANGES
Rev.C	06/2023	New Issue for IVDR Regulation (UE) 2017/746 compliance