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
Intended use

In vitro diagnostic medical device for performing test-tube agglutination tests for titration of anti-brucella antibodies according to Wright. For professional use only.

Principle of the Method

When serum containing specific agglutinins reacts with homologous antigen under optimized controlled conditions it is able to cause a visible agglutination. The degree of agglutination depends on the concentration of the antigen, antibodies, salts in the medium and on the temperature.

Storage and stability

Storage temperature 2-8 °C

Stored at 2-8°C, avoiding direct light, the reagents are stable until the expiry date stated on the label. Do not freeze. Avoid microbial contamination. Repeated stability tests on three different batches of reagent confirmed the validity of the reagent for at least 36 months at 2-8°C. A slight variation in the composition of the reagents can occur from batch to batch, without interfering with the test results

PRECAUTIONS and WARNINGS

- Reagents and waste materials shall be disposed of in accordance with Community waste provisions or national or regional provisions.
- Reagents may contain non-active components such as preservatives and detergents. The total concentration of these components is below the limits set out in Regulation 1272/2008 EC and subsequent amendments and additions.
- It is recommended that the reagent be handled in accordance with the rules of good laboratory practice and that appropriate personal protective equipment be used.

- Do not use the reagent if it is visibly degraded (e.g. presence of corpuscles).
- All human samples shall be handled and disposed of as potentially infectious material.
- The kit should only be used by qualified and properly trained technical personnel.
- Diagnoses shall be carried out exclusively by authorised and qualified personnel.
- Comply with national directives on occupational safety and quality assurance.
- Use equipment that complies with current regulations.

Reporting of serious incidents

In the event of a serious incident in relation to the use of the device, please inform the manufacturer (via your distributor) and the competent authority of the European Union member state where the incident occurred. For other jurisdictions, reporting must be made in accordance with regulatory requirements. Reporting serious incidents helps provide more information about the safety of the diagnostic medical device

PROCEDURE

Quality control Kit suspensions should be analyzed with positive and negative control serum. The absence of positive or negative reactions respectively indicates alteration of suspensions and/or controls.

Sample: Clear serum, diluted 1:25 with Saline Solution

Diagnostic suspension: Reagents and materials must be brought to room temperature before use.

Shake the antigen bottles thoroughly until a homogeneous suspension is obtained.

Tube	Nr	1	2	3	4	5	6	K+	K-
Saline	mL	-	0,5	0,5	0,5	0,5	0,5	-	-
Diluted serum 1:25	mL	0,5	0,5	→	→	→		-	-
Positive Control	mL	-	-	-	-	-	-	0,5	-
Negative control	mL	-	-	-	-	-	-	-	0,5
Perform doubling dilutions by transferring 0,5 mL from Test Tube N.2 into the next, to Test tube N.6, discarding the last 0,5 mL remaining.									
Diagnostic suspension	mL	0,5	0,5	0,5	0,5	0,5		0,5	0,5
Shake and incubate in a thermostat at 37°C for 24 hours									
Titer		50	100	200	400	800	1600	-	-

Reading: Without shaking the tubes, the possible degree of sedimentation and the clarification of the supernatant is observed. The control tube must not show evident agglutination.

The titre of antibodies is given by the reciprocal of the highest dilution of the sample showing a visible agglutination. Titres between 50 and 100 must be considered as suspect; titres of 100 or higher are confirmatory for diagnosis of disease.

Validation test:

Specificity: Tests performed with three different lots of Brucella bacterial suspensions on negative human samples had repeatedly shown negative results.

Precision: Repeatability tests (in the run) and Reproducibility tests (between run) performed with three different lots of Brucella bacterial suspensions, on human samples and animal sera with known titre, had constantly shown the expected results ± 1 dilution (Method accuracy limit)

REVISIONE/ REVISION	DATA/DATE	MODIFICHE/CHANGES
C	05/2023	Nuova emissione per adeguamento IVDR Regolamento (UE) 2017/746/ New issue for IVDR Regulation (UE) 2017/746 compliance

- Do not use the reagent if it is visibly degraded (e.g. presence of corpuscles).
- All human samples shall be handled and disposed of as potentially infectious material.
- The kit should only be used by qualified and properly trained technical personnel.
- Diagnoses shall be carried out exclusively by authorised and qualified personnel.
- Comply with national directives on occupational safety and quality assurance.
- Use equipment that complies with current regulations.

Reporting of serious incidents

In the event of a serious incident in relation to the use of the device, please inform the manufacturer (via your distributor) and the competent authority of the European Union member state where the incident occurred. For other jurisdictions, reporting must be made in accordance with regulatory requirements. Reporting serious incidents helps provide more information about the safety of the diagnostic medical device

PROCEDURE

Quality control Kit suspensions should be analyzed with positive and negative control serum. The absence of positive or negative reactions respectively indicates alteration of suspensions and/or controls.

Sample: Clear serum, diluted 1:25 with Saline Solution

Diagnostic suspension: Reagents and materials must be brought to room temperature before use.

Shake the antigen bottles thoroughly until a homogeneous suspension is obtained.

Tube	Nr	1	2	3	4	5	6	K+	K-
Saline	mL	-	0,5	0,5	0,5	0,5	0,5	-	-
Diluted serum 1:25	mL	0,5	0,5	→	→	→		-	-
Positive Control	mL	-	-	-	-	-	-	0,5	-
Negative control	mL	-	-	-	-	-	-	-	0,5
Perform doubling dilutions by transferring 0,5 mL from Test Tube N.2 into the next, to Test tube N.6, discarding the last 0,5 mL remaining.									
Diagnostic suspension	mL	0,5	0,5	0,5	0,5	0,5		0,5	0,5
Shake and incubate in a thermostat at 37°C for 24 hours									
Titer		50	100	200	400	800	1600	-	-

Reading: Without shaking the tubes, the possible degree of sedimentation and the clarification of the supernatant is observed. The control tube must not show evident agglutination.

The titre of antibodies is given by the reciprocal of the highest dilution of the sample showing a visible agglutination. Titres between 50 and 100 must be considered as suspect; titres of 100 or higher are confirmatory for diagnosis of disease.

Validation test:

Specificity: Tests performed with three different lots of Brucella bacterial suspensions on negative human samples had repeatedly shown negative results.

Precision: Repeatability tests (in the run) and Reproducibility tests (between run) performed with three different lots of Brucella bacterial suspensions, on human samples and animal sera with known titre, had constantly shown the expected results ± 1 dilution (Method accuracy limit)

REVISIONE/ REVISION	DATA/DATE	MODIFICHE/CHANGES
C	05/2023	Nuova emissione per adeguamento IVDR Regolamento (UE) 2017/746/ New issue for IVDR Regulation (UE) 2017/746 compliance