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Formato	REF	Sospensione	Composizione
Kit 5 x 100 test	CSI087310	Antigen Set	5 flaconi x 5 ml
Kit 3 x 100 test	CSI087230	Brucella	3 flaconi x 5 ml
Kit 3 x 100 test	CSI087235	S. typhi (totale)	3 flaconi x 5 ml
Kit 3 x 100 test	CSI087234	S. typhi (antigene O)	3 flaconi x 5 ml
Kit 3 x 100 test	CSI087233	S. typhi (antigene H)	3 flaconi x 5 ml
Kit 3 x 100 test	CSI087232	S. paratyphi A (totale)	3 flaconi x 5 ml
Kit 3 x 100 test	CSI087311	S. paratyphi A (antigene O)	3 flaconi x 5 ml
Kit 3 x 100 test	CSI087312	S. paratyphi A (antigene H)	3 flaconi x 5 ml
Kit 3 x 100 test	CSI087231	S. paratyphi B (totale)	3 flaconi x 5 ml
Kit 3 x 100 test	CSI087313	S. paratyphi B (antigene O)	3 flaconi x 5 ml
Kit 3 x 100 test	CSI087314	S. paratyphi B (antigene H)	3 flaconi x 5 ml
Kit 3 x 100 test	CSI087236	S. paratyphi C (totale)	3 flaconi x 5 ml
Kit 3 x 100 test	CSI087237	S. paratyphi C (antigene O)	3 flaconi x 5 ml
Kit 3 x 100 test	CSI087238	S. paratyphi C (antigene H)	3 flaconi x 5 ml
-	CSI087388	Febbrile Controllo Negativo	1 flacone x 2 ml

DESTINAZIONE D'USO

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Bibliografia

- Robert W. Dornier et al. Clinica Chimica Acta 1987; 167: 1-21
- Frederick Wolfe et al. Arthritis and Rheumatism 1991; 34: 951-960
- Robert H Shmerling et al. The American Journal of Medicine 1991; 91: 528-534
- Koriz T N et al. Journal of Immunological Methods 1980; 32: 1-9
- Assameh S N et al. Journal of Immunological Methods 1980; 34: 205-215
- Young DS Effects of drugs on clinical laboratory test, 4th ed. AACCC Press, 1995.

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INTENDED USE

Colored bacterial suspensions for titration on microplate of human serum antibodies. The results of the test must always be interpreted in conjunction with the clinical context. For professional use only.

CLINICAL SIGNIFICANCE

Serological diagnosis of infectious diseases characterized by persistent fever is dependent upon demonstration of agglutination reaction between the appropriate antigen and the patient's serum containing antibodies. Grumbaum and Widal first introduced routine immunological applications. Their test has become universally known as the "Widal test" and quantitatively measures the antibodies (agglutinins) in the sera of patients with Typhoid Fever. Furthermore, the use of serological agglutination tests has become quite common in assisting in the diagnosis of brucellosis in man and other animals.

PRINCIPLE OF THE METHOD

When serum containing specific agglutinins reacts with homologous antigen under optimized conditions, it is able to cause a visible agglutination. This reaction may be performed in test tubes, on slides or in microplate. In Plate test, the two-fold diluted serum is dispensed in microplate wells and after antigen addition, the reaction is incubated for the stated time and read. The microorganisms are dyed to facilitate the reading. Positive reactions appear as a smooth mat of cells on the bottom of wells, while negative reactions show a totally or partially compact button of colored cells.

COMPONENTS

All concentrations refer to ready to use reagent:

Suspension (the specific name is reported on vial label)
Dyed bacterial suspension at optimal concentration for microplate test.
Contains Sodium azide 0,95 g/L.

Positive Control Serum (polyvalent, animal source)
Contains Sodium azide 0,95 g/L.

***Warning:** Products that contain sodium azide may react with lead and copper in plumbing to form explosive metal azides. On disposal flush with large volumes of water to prevent azide build up.

4 96-well microplates U bottom wells to perform analysis.

Other material required, but not supplied

Physiological saline, automatic pipettes to distribute 10 µL, 50 µL, 90 µL.

PRECAUTIONS and WARNINGS

- Reagents and waste materials shall be disposed of in accordance with Community waste provisions or national or regional provisions.
- Reagents may contain non-active components such as preservatives and detergents. The total concentration of these components is below the limits set out in Regulation 1272/2008 EC and subsequent amendments and additions.
- It is recommended that the reagent be handled in accordance with the rules of good laboratory practice and that appropriate personal protective equipment be used.
- Do not use the reagent if it is visibly degraded (e.g. presence of corpuscles).
- All human samples shall be handled and disposed of as potentially infectious material.
- The kit should only be used by qualified and properly trained technical personnel.
- Diagnoses shall be carried out exclusively by authorised and qualified personnel.
- Comply with national directives on occupational safety and quality assurance.
- Use equipment that complies with current regulations.

Reporting of serious incidents

In the event of a serious incident in relation to the use of the device, please inform the manufacturer (via your distributor) and the competent authority of the European Union member state where the incident occurred. For other jurisdictions, reporting must be made in accordance with regulatory requirements. Reporting serious incidents helps provide more information about the safety of the diagnostic medical device.

REAGENT PREPARATION

Reagents are liquid ready to use. Avoid any chemical or bacteriological contamination. Mix well reagents before use. Suspension must appear uniform and free of visible particles.

Preservation and Stability

 = Storage Temperature 2-8 °C

Store reagents and controls at 2-8° C avoiding direct light. Do not freeze. Under above conditions reagents are stable until expiration date indicated on the label. Slight variations in composition among batches will not affect test results.

SPECIMEN COLLECTION AND PRESERVATION

Not-inactivated sera must be used. Samples can be stored at 2-8° C for 48 h before testing. For longer time of storage samples must be frozen (-20° C). Highly hemolyzed, lipemic or contaminated samples cannot be used. Samples with the presence of fibrin should be centrifuged before testing.

PROCEDURE

Quality Control

Kit suspensions should be tested with positive control serum and with saline as negative control.

Failure to obtain a positive and negative result with the respective controls evidence a deterioration of suspension and/or controls.

Technique

- Bring reagents and materials to room temperature. Low temperature can reduce the method sensitivity.
- Swirl gently the reagent before using. Suspension must appear uniform and free from visible particles.

		ROWS											
		1A	2A	3A	4A	5A	6A	7A	8A	9A	10A	11A	12A
COLUMNS	1B	2B	3B	4B	5B	6B	7B	8B	9B	10B	11B	12B	
	1C	2C	3C	4C	5C	6C	7C	8C	9C	10C	11C	12C	
	Etc	-	-	-	-	-	-	-	-	-	-	-	-

3 - Starting from the first well of the first column (1 A), dilute 10 µL of sample with 90 µL saline.

4 - Distribute 50 µL of saline in all the wells of the row, until the 10 A; in 11 A add 50 µL of the Positive control and in 12 A, 50 µL of negative control.

5 - Dilute the sample transferring 50 µL from the first to the second well.

6 - Mix by pipetting several times and transfer 50 µL to well 3A.

7 - Follow up to well 10 A discarding the last 50 µL.

8 - Repeat this operation for each sample.

9 - Distribute 50 µL of bacterial suspension in each well.

10 - Shake for 30 seconds with a plate mixer or manually.

11 - Incubate at 37° C in a moist chamber in absence of vibrations for 16 ± 4 hours.

To improve the reading (especially for Brucella) a further incubation of 30 minutes at room temperature on a surface without vibrations is suggested.

Results

Positive Reaction (agglutination): Homogeneous layer of cells without sedimented bacteria (clear supernatant).

Negative reaction (absence of agglutination): sedimented cells forming a button on the well bottom.

Sometimes with positive control it is possible to observe traces of agglutinated bacteria on the bottom of the well. These are not to be confused with the button typical of the negative reaction.

Titer

The titer is defined as reciprocal dilution of serum showing a visible agglutination.

Results interpretation

Generally, titers in the range 80-160 are considered suspect; only higher titers are probative for the diagnosis of the disease. Significant titers may be obtained from individuals immunized with typhoid vaccine. Must be noted that immunological response for each subject against an infectious bacterial agent could be influenced by several factors.

Warnings

Prozone reactions, especially with Brucella reagent, are possible. If clumping is seen in lower dilutions, care should be taken to observe agglutination at the higher dilutions (at least until 1:640). Serum agglutination apparently negative may be due to an excess of antibodies, which, at low dilutions, prevents the formation agglutinate (prozone phenomenon). In case of not significant positivity, the addition of serum anti - human Ig may reveal the presence of incomplete antibodies (Coombs test) present in particular in the case of chronic disease. As with all serological procedures, a single antibody determination should not be used for diagnosis. The significance of antibody presence, in a single serum can provide only the evidence of prior exposure to etiologic agent. To give a correct diagnosis it is necessary to observe a significant increase of titer between samples collected 10-14 days apart.

Method features

Results obtained by Plate test showed a good correlation (±1 dilution) with titers obtained by classical tube tests.

Validation Test

Specificity





Tests carried out with three different lots of bacterial suspension "O", "H", "O+H", Brucella on human negative sera had given repeatedly negative results.

Sensitivity

Tests carried out with three different lots of bacterial suspension "O", "H", "O+H", Brucella on human positive sera, had confirmed the titer obtained with reference methods.

Stability

Stability test, in real time, was carried out on three different lots of each bacterial suspension have confirmed the functionality of the reagent for at least 36 months to 2-8°C.

Symbols used in IFU and Packaging	
 In vitro diagnostic medical device vitro	 Manufacturer
 Catalogue Number	 Instruction for use
 Lot Number	 Temperature limitation
 Expiration date	

Bibliography

- Robert W. Dörner et al. Clinica Chimica Acta 1987; 167: 1-21
- Frederick Wolfe et al. Arthritis and Rheumatism 1991; 34: 951-960
- Robert H Shmerling et al. The American Journal of Medicine 1991; 91: 528-534
- Koriz T N et al. Journal of Immunological Methods 1980; 32: 1-9
- Assameh S N et al. Journal of Immunological Methods 1980; 34: 205-215
- Young DS Effects of drugs on clinical laboratory test, 4th ed. AACCC Press, 1995.

REVISION	DATE	CHANGE
F	04/2023	New Issue for IVDR Regulation (UE) 2017/746 compliance