LOT	Lot: 25/001		Exp: 31/08/2025		Storage: 2-8 °C
------------	--------------------	--	------------------------	--	------------------------

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Bibliografia

Manuale "Biosafety in Microbiological and biomedical Laboratories", Centers for Disease Control / National Institute of Health, 1984.

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INTENDED USE

the serum "ISOCONTROL Serum CK" is used as a quality control to verify the accuracy and precision of the method quantitative immunochemical diagnostic Sclavo. The control values are described in the theoretical section "Assigned values". For professional use only.

REAGENT INFORMATION

Reagent: n° 2 vial x 5 mL

The serum "ISOCONTROL Serum CK" is obtained from defibrinated human plasma and freeze-dried with the addition of recombinant CK. The nature of these isoenzymes allows a great degree of affinity with the sera of the patients concerned.

Reconstitution of serum

Remove the cap carefully and avoid wasting traces of lyophilized. Add 5 mL of distilled water exactly, cap, invert a few times (do not shake) to avoid foaming, leave to stand for about 30 minutes.

Assigned values

Analyte	Assigned Value		
	Value	Range	U.M.
CK NAC	220	165-275	U/L
CK MB	52	39-65	U/L

Storage and stability

= Storage temperature 2-8°C

If stored closed at 2-8 °C, avoiding direct light, are stable until the expiration date, printed on the label. Stability after reconstitution: 4 hours at 18-25°C; 7 days at 2-8 °C (With closed vial); 30 days at -15-25°C

PRECAUTIONS and WARNINGS

1. Reagents and waste materials shall be disposed of in accordance with Community waste provisions or national or regional provisions.
2. In addition to any risk claims relating to active components, reagents may contain non-active components such as preservatives and detergents. The total concentration of these components is below the limits set out in Regulation 1272/2008 EC and subsequent amendments and additions.
3. It is recommended that the reagent be handled in accordance with the rules of good laboratory practice and that appropriate personal protective equipment be used.
4. Do not use the reagent if it is visibly degraded (e.g. presence of corpuscles) After reconstitution
5. All human samples shall be handled and disposed of as potentially infectious material. The sera have been tested and proved negative for the presence of HbsAg, and also for anti-HIV 1 and 2 and anti-HCV antibodies. As there is no possibility to exclude definitely that products derived from human blood transmit infectious agents, it is recommended to handle the control with the same precautions used for patient specimens.
6. The kit should only be used by qualified and properly trained technical personnel.
7. Diagnoses shall be carried out exclusively by authorised and qualified personnel.
8. Comply with national directives on occupational safety and quality assurance.
9. Use equipment that complies with current regulations.

Reporting of serious incidents

Please inform the manufacturer (through your distributor) and the competent authority of the member state of the European Union in which the user and/or patient is established, of cases of serious incident that has occurred in relation to the device. For other jurisdictions, reports of serious incidents must be made in accordance with the regulatory requirements of the home Member State. By reporting serious incidents, you help provide more information about the safety of your in vitro medical diagnostic device.

Values assignment

The theoretical values were calculated on the basis of the results on automatic biochemical analyzers and reagents Sclavo Diagnostics Sclavo Diagnostics.

Materials not included in the kit

- Reagents for analytical systems Sclavo Diagnostics

Symbols used in IFU and Packaging	
In vitro Diagnostics Medical Device	Manufacturer
Catalogue Number	Instructions for Use
Lot Number	Storage Temperature
Expiration Date	Biological risk

Bibliography

Manual "Biosafety in Microbiological and biomedical Laboratories", Centers for Disease Control/National Institute of Health, 1984.

REVISION	DATE	CHANGE
Rev. E	11/2024	Packaging Format change; reconstitution of serum and Storage and stability