|                               |                                |                               |
|-------------------------------|--------------------------------|-------------------------------|
| <b>LOT</b> Lot: <b>24/001</b> | <b>Exp.:</b> <b>21/11/2025</b> | <b>Storage:</b> <b>2-8 °C</b> |
|-------------------------------|--------------------------------|-------------------------------|

#### ORDER INFORMATION

| Format       | Code                 | Composition       |
|--------------|----------------------|-------------------|
| Kit 3 x 1 mL | <b>REF</b> B35282591 | n° 3 vials x 1 mL |

#### Intended use

To perform the calibration for the determination of total protein in human urine, in association with the kit and the total urinary protein determined by instrumental applications Sclavo Diagnostics International.

#### Features of the calibrator

The calibrator consists of human proteins, added in defined amounts to a protein solution buffered and stabilized.

#### Composition

The calibration set consists of three calibrators with the exact concentration of protein. The analyte concentration was determined by quantitative methods using reagents Sclavo Diagnostics Int.

The concentrations are lot dependent; the exact values are specified in the table below. The concentrations are expressed in "mg/dL".

| Analyte | Assigned Value | U.M   |
|---------|----------------|-------|
| PTU     | <b>91.0</b>    | mg/dL |

#### Storage and stability of the serum

The calibrators are liquid, ready to use. Are stable until the expiration date printed on the individual bottles if stored at a temperature between 2 and 8 °C

#### PRECAUTIONS AND WARNINGS

1. Reagents and waste materials shall be disposed of in accordance with Community waste provisions or national or regional provisions.
2. In addition to any risk claims relating to active components, reagents may contain non-active components such as preservatives and detergents. The total concentration of these components is below the limits set out in Regulation 1272/2008 EC and subsequent amendments and additions.
3. It is recommended that the reagent be handled in accordance with the rules of good laboratory practice and that appropriate personal protective equipment be used.
4. Do not use the reagent if it is visibly degraded (e.g. presence of corpuscles) after reconstitution.
5. The human source material used for this product has been tested and found negative (or not repeatedly reactive) for the presence of antibodies to HIV and HCV as well as for hepatitis B surface antigen. Since no known method can ensure that infectious agents are absent, these reagents should be handled with the precautions required for each sample of serum or potentially infected blood.

#### Reporting of serious incidents

Please inform the manufacturer (through your distributor) and the competent authority of the member state of the European Union in which the user and/or patient is established, of cases of serious incident that has occurred in relation to the device. For other jurisdictions, reports of serious incidents must be made in accordance with the regulatory requirements of the home Member State. By reporting serious incidents, you help provide more information about the safety of your in vitro medical diagnostic device.

#### Assigned values and traceability

The theoretical values have been calculated on the basis of the results obtained on Sclavo Diagnostics biochemical analyzers using Sclavo Diagnostics reagents. The traceability of the values, with reference to each diagnostic method, is reported on the enclosed product information sheet. The use of procedures or reagents other than those indicated may lead to differences between the expected values and those obtained.

#### Materials required but not included in the kit

- Distilled or deionized water
- Reagents for Sclavo Diagnostics analytical systems

| Symbols used in IFU and Packaging                                                                   |                                                                                                         |
|-----------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------|
| <b>IVD</b> In vitro diagnostic medical device vitro                                                 |  Manufacturer        |
| <b>REF</b> Catalogue Number                                                                         |  Instruction for use |
| <b>LOT</b> Lot Number                                                                               |  Storage Temperature |
|  Expiration date |  Biological Risk     |

#### Reference

Data about the village of File Sclavo Diagnostics International

| REVISION | DATE    | CHANGE                                                 |
|----------|---------|--------------------------------------------------------|
| Rev.D    | 06/2023 | New Issue for IVDR Regulation (UE) 2017/746 compliance |

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