CONTROLLO LUE / LUE CONTROL

CSI053060

1x1 ml

**Intended use**

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Declaration of Serious incident

Please notify to the manufacturer (through your distributor) and the competent authority of the Member state of the European union in which the user or patient are established, any serious incident that as occurred in relation of this device. By reporting serious incidents, you provide information that can contribute to the safety of in vitro medical devices

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RPR:

Symbols used in IFU and Packaging

In vitro diagnostic medical device vitro	Manufacturer
Catalogue Number	Instruction for use
Lot Number	Temperature limitation
Expiration date	

REVISION	DATE	CHANGE
Rev.B	08/2024	New format

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1x1 ml

**Intended use**

Freeze dried Human positive serum to be used as Quality Control in the Syphilis IVD Tests

Method of use

The Lue Control is freeze dried. Before use, remove the cap carefully, to avoid the loss of part of the lyophilized serum. Add exactly 1 mL of distilled water, plug, mix by inversion and let it stand for about 30 minutes. Again, swirl gently to avoid foam production. Avoid bacterial contaminations. On the vial label is reported the expiring date of the unopened vial.

Warning: Biohazard

Each donor unit, used in the preparation of this material, has been tested and found to be negative for the presence of the antibodies to HIV and HCV as well as for hepatitis B surface antigen. As it's impossible to guarantee the lack of infectious agents, these materials should be handled, as recommended, as for any potentially infectious human material Biosafety Level 2)*.

* For more informations is we suggest to look the handbook "Biosafety in Microbiological and Biomedical Laboratories", Centers for Disease Control/National Institutes of Health, 1984.

Ingredients

Human serum positive for anti-Syphilis antibodies.

Storage

Store at 2-8 ° C avoiding exposure to light. Do not use if the product appears to be deteriorated.

Stability

The stability of the product is printed on the bottle.

Characteristics and stability of the reconstituted serum are the same of the fresh serum. If required from the method, before use inactivate the serum at 56°C for 30 min.

Declaration of Serious incident

Please notify to the manufacturer (through your distributor) and the competent authority of the Member state of the European union in which the user or patient are established, any serious incident that has occurred in relation of this device. By reporting serious incidents, you provide information that can contribute to the safety of in vitro medical devices

Disposal

Disposal must be performed in accordance with the EC regulations regarding waste, or the local national or regional legislation.

Titre:

V.D.R.L.

Complement fixation (according to Kolmer)

- Cardiolipinic Antigen:
- Treponemic Antigen:

TPHA:

RPR:

Symbols used in IFU and Packaging	
In vitro diagnostic medical device vitro	Manufacturer
Catalogue Number	Instruction for use
Lot Number	Temperature limitation
Expiration date	

REVISION	DATE	CHANGE
Rev.B	08/2024	New format

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Sclavo Diagnostics International S.p.A : - Loc.Pian dei Mori, Via Po 26-28 - 53018 Sovicille (Si) Italy
Sclavo Diagnostics International S.p.A. - Loc.Pian dei Mori, Via Po 26-28 - 53018 Sovicille (Si) Italy