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### INTENDED USE

The detergent (Methyl Alcohol) can be used manually in the preparation of the slide in the fixation phase or on automatic systems.

### METHOD OF USE

The manual application, in the slide preparation and fixing phases, is described in the technical data sheet supplied with the stains. The method of use on automatic systems is described in the relative operating manuals.

### INFORMATION ON REAGENT

Sysclean		
Methyl Alcohol	95%	
Acetone	5%	

#### \*Warning: DANGER

Contains: Methanol (CAS 67-56-1)

**H225** - Highly flammable liquid and vapor.

**H301+H311+H331** - Toxic if swallowed, in contact with skin or if inhaled

**H370** - Causes damage to organs.

**P210** - Keep away from heat sources, hot surfaces, sparks, open flames or other sources of ignition. Not smoking.

**P280** - Wear protective gloves/ protective clothing/eye protection/face protection/ hearing protection/...

**P308+P311** - IF exposed or concerned: Call a POISON CENTER or a doctor

**P403+P233** - Store in a well-ventilated place. Keep container tightly closed.

**P370+P378**- In case of fire: extinguish with suitable extinguishing means (CO2/POWDER). Do not use water

### Storage

The reagent is stable up to the expiry date indicated on the tank if stored at a temperature between 15 and 25 ° C

### Stability

The pad, after the first opening, must be closed and stored at a temperature between +15°C and + 25 °C, under these conditions it is declared stable until the expiration date described on the external label of the container.

### PRECAUTIONS and WARNINGS

1. The kit must be used only by qualified and adequately trained technical personnel.
2. Diagnoses should only be carried out by authorized and qualified personnel.
3. It is recommended to handle the reagent according to good laboratory practice and to use suitable personal protective equipment.
4. Observe the national regulations on occupational safety and quality assurance.
5. Laboratory standards relating to infection protection should be used

### Instructions for disposal

The used solutions and the expired solutions must be disposed of as hazardous waste in accordance with the Community provisions on waste or with the national or regional provisions in force.

### Classification of dangerous substances

The classification, labelling and packaging of products comply with EC regulation 1272/2008 and subsequent amendments and additions. According to this regulation, the reagent is classified as DANGEROUS and must be used with caution. Observe the classification of hazardous substances on the label and follow the instructions in the Safety Data Sheet. The safety data sheet is available on the website or on request.

### Reporting of serious incidents

Please inform the manufacturer (through its distributor) and the competent authority of the member state of the European Union in which the user and / or patient is established, of the cases of serious accident occurring in relation to the device. For other jurisdictions, reporting of serious incidents must be made in accordance with the regulatory requirements of the home member state. By reporting serious incidents, you help provide more information about the safety of the in vitro diagnostic medical device.

Symbols used in IFU and Packaging	
 In vitro diagnostic medical device vitro	 Manufacturer
 Catalogue Number	 Instruction for use
 Lot Number	 Temperature limitation
 Expiration date	

REVISION	DATE	CHANGE
Rev.C	01/2025	Update of the "Method of use" chapter.