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INTENDED USE

The Buffer can be used with the manual method in slide preparation during the staining phase, or in automatic systems in the staining phase. For professional use only.

METHOD OF USE

The manual application, in the slide preparation and staining phases, is described in the technical data sheet supplied with the stains. The automatic method on "Sysmex SP1000 and SP100" systems is described in the relative operating manuals.

PRINCIPLE OF THE METHOD

pH 6.8 buffers are intended to create and maintain an environment with a constant pH of 6.8 in vitro, ensuring that chemical and biological reactions take place under controlled and predictable conditions. This buffer solution is a fundamental tool for the research and analysis of samples subject to histological staining.

INFORMATION ON REAGENT

Buffer Solution	
Phosphate Buffer	pH 6.8

*Warning: ATTENTION

Contains: 5-chloro-2-methyl-2H-isothiazol-3-one (EC 247-500-7); 2-methyl-2H-isothiazol-3-one (EC 220-239-6); (3:1) May cause an allergic reaction.

Conservation

The reagent is stable until the expiry date printed on the tank's label when stored at a temperature between 15 and 25 ° C

Stability

The buffer, after the first opening, must be closed and stored at a temperature between +15°C and +25°C. Under these conditions it is declared stable until the expiry date described on the external label of the container.

Instructions for disposal

The used solutions and the expired solutions must be disposed of as hazardous waste in accordance with the Community provisions on waste or with the national or regional provisions in force.

Reporting of serious incidents

In the event of a serious incident in relation to the use of the device, please inform the manufacturer (via your distributor) and the competent authority of the European Union member state where the incident occurred. For other jurisdictions, reporting must be made in accordance with regulatory requirements. Reporting serious incidents helps provide more information about the safety of the diagnostic medical device.

PRECAUTIONS and WARNINGS

- The kit should only be used by qualified and properly trained technical personnel.
- Diagnoses shall only be carried out by authorised, qualified and appropriately trained personnel.
- It is recommended to handle the reagent in accordance with the rules of good laboratory practice and to use appropriate personal protective equipment.
- Comply with national directives on occupational safety and quality assurance.
- Laboratory standards for protection against infection shall be used.
- Do not dispose of the product with household waste Do not enter into the sewers.
- Reagents and waste materials shall be disposed of in accordance with Community waste provisions or national or regional provisions.
- Reagents may contain non-active components such as preservatives and detergents. The total concentration of these components is lower than the limits set out in Regulation 1272/2008 EC and subsequent amendments and additions.
- Do not use the reagent if it is visibly degraded (e.g. presence of corpuscles).
- Comply with national directives on occupational safety and quality assurance.
- Use equipment that complies with current regulations.

Symbols used in IFU and Packaging	
 In vitro diagnostic medical device vitro	 Manufacturer
 Catalogue Number	 Instruction for use
 Lot Number	 Temperature limitation
 Expiration date	

REVISION	DATE	CHANGE
Rev. A	10/2023	First issue