<b>LOT</b> Lot: <b>23/001</b>	<b>Exp.:</b> <b>2026/08/31</b>	<b>Conservazione / Storage: 2-8 °C</b>
-------------------------------	--------------------------------	--

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#### RECOMMENDED USE

Product for use in the quality control and determination of the concentration of HDL and LDL in human serum or plasma.

#### REAGENT INFORMATION

**Reagent:** n° 5 vials x 3 mL freeze-dried Code 191153

Control serum is obtained from defibrinated freeze-dried human plasma. The declared nature of the product denotes the potential risk involved and it must therefore be used adopting same precautions used in handling test samples.

#### Instructions for use

- Resuspend the vial content with con 1 ml of Distilled water
- Replace the cap and mix gently by inversion, avoiding foam formation
- Mix gently and let stand before use: 30 minutes for HDL Cholesterol and 2 hours for LDL

#### Assigned value

Analyte	Assigned value		Unit
	Value	Range	
Cholesterol HDL	<b>54</b>	<b>41 - 68</b>	<b>mg/dL</b>
	1.40	1.12 - 1.68	mmol/L
Cholesterol LDL	<b>105</b>	<b>78 - 131</b>	<b>mg/dL</b>
	2.70	2.30 - 3.11	mmol/L

#### Conservation

Unopened vials are stable until the expiration date printed on the individual bottles if stored at a temperature between 2 and 8 °C.

Reconstituted open vials are stable:

8 hours at 25°C - 7 days at 2-8°C - 30 days at -20°C

#### PRECAUTIONS and WARNINGS

1. Reagents and waste materials shall be disposed of in accordance with Community waste provisions or national or regional provisions.
2. In addition to any risk claims relating to active components, reagents may contain non-active components such as preservatives and detergents. The total concentration of these components is below the limits set out in Regulation 1272/2008 EC and subsequent amendments and additions.
3. It is recommended that the reagent be handled in accordance with the rules of good laboratory practice and that appropriate personal protective equipment be used.
4. Do not use the reagent if it is visibly degraded (e.g. presence of corpuscles) After reconstitution
5. All human samples shall be handled and disposed of as potentially infectious material. The sera have been tested and proved negative for the presence of HbsAg, and also for anti-HIV 1 and 2 and anti-HCV antibodies. However, the material must be considered potentially dangerous and handled in the same way as the test samples.
6. The kit should only be used by qualified and properly trained technical personnel.
7. Diagnoses shall be carried out exclusively by authorised and qualified personnel.
8. Comply with national directives on occupational safety and quality assurance.
9. Use equipment that complies with current regulations.

#### Reporting of serious incidents

Please inform the manufacturer (through your distributor) and the competent authority of the member state of the European Union in which the user and/or patient is established, of cases of serious incident that has occurred in relation to the device. For other jurisdictions, reports of serious incidents must be made in accordance with the regulatory requirements of the home Member State. By reporting serious incidents, you help provide more information about the safety of your in vitro medical diagnostic device.

#### Values assignment

The theoretical values were calculated on the basis of the results on automatic biochemical analyzers and reagents Sclavo Diagnostics Sclavo Diagnostics.

#### Materials not included in the kit

- Reagents for analytical systems Sclavo Diagnostics

Symbols used in IFU and Packaging	
In vitro diagnostic medical device vitro	Manufacturer
Catalogue Number	Instruction for use
Lot Number	Temperature limitation
Expiration date	Biological Risk

REVISION	DATE	CHANGE
Rev. F	01/2023	Revision for compliance with IVDR Regulation (UE) 2017/746