<b>LOT</b> Lot:	<b>25/001</b>	 Exp.:	<b>2026/12/26</b>	 Storage:	<b>2-8 °C</b>
-----------------	---------------	---	-------------------	---	---------------

#### Intended Use

Kit to perform the accuracy control for determination of Total Protein in human urine, in association with the kit Total Urinary Protein SDS determined by instrumental applications Scclavo Diagnostics International.

#### Characteristics of controls

The controls consist of human proteins, added in defined amounts to a buffered and stabilized solution.

#### Kit Composition

The set of controls consists of three High control with exact concentration of protein and Low in three controls the exact concentration of protein.

The concentration of the analyte was determined by quantitative methods using reagents Scclavo Diagnostics Int. The concentrations are lot dependent; the exact values are specified in the table below and are expressed in "mg/dL".

Controllo LOW			
Analyte	Assigned Value	Range	U.M.
PTU	23	17-28	mg/dL
Controllo HIGH			
Analyte	Assigned Value	Range	U.M.
PTU	167	125-209	mg/dL

#### Storage and stability of the serum

The controls are liquid, ready to use.

Are stable until the expiration date printed on the individual bottles if stored at a temperature between 2 and 8 ° C in the original package.

#### PRECAUTIONS and WARNINGS

1. Reagents and waste materials shall be disposed of in accordance with Community waste provisions or national or regional provisions.
2. In addition to any risk claims relating to active components, reagents may contain non-active components such as preservatives and detergents. The total concentration of these components is below the limits set out in Regulation 1272/2008 EC and subsequent amendments and additions.
3. It is recommended that the reagent be handled in accordance with the rules of good laboratory practice and that appropriate personal protective equipment be used.
4. Do not use the reagent if it is visibly degraded (e.g. presence of corpuscles) After reconstitution
5. Caution – Biohazard. The human source material used for this product has been tested and found negative (or not repeatedly reactive) for the presence of antibodies to HIV and HCV as well as for hepatitis B surface antigen. Since no known method can ensure that infectious agents are absent, these reagents should be handled with the precautions required for each sample of serum or potentially infected blood.

#### Reporting of serious incidents






Please inform the manufacturer (through your distributor) and the competent authority of the member state of the European Union in which the user and/or patient is established, of cases of serious incident that has occurred in relation to the device. For other jurisdictions, reports of serious incidents must be made in accordance with the regulatory requirements of the home Member State. By reporting serious incidents, you help provide more information about the safety of your in vitro medical diagnostic device.

#### Assigned values and traceability

The theoretical values have been calculated on the basis of the results obtained on Scclavo Diagnostics biochemical analyzers using Scclavo Diagnostics reagents. The use of procedures or reagents other than those indicated may lead to differences between the expected values and those obtained.

#### Materials required but not included in the kit

- Reagents for Scclavo Diagnostics analytical systems

Symbols used in IFU and Packaging	
 In vitro diagnostic medical device vitro	 Manufacturer
 Catalogue Number	 Instruction for use
 Lot Number	 Temperature limitation
 Expiration date	 Biological Risk

REVISION	DATE	CHANGE
Rev. D	03/2024	Revision for compliance with IVDR Regulation (UE) 2017/746

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