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INTRODUCTION

Gynecult broth is intended for the isolation and cultivation of microorganisms of the genus *Trichomonas* and yeasts from urogenital samples. *Trichomonas vaginalis* and yeasts are the most common causative agents of vaginitis (25-50% of infected women has symptoms as dysuria, itch and vaginal heartburn). *T. vaginalis* can be also present in the male urogenital system but in this case the infection is asymptomatic.

PRINCIPLE OF THE METHOD

The broth contains pepton, cystein, yeast extract as aminoacid, nitrogen sulfur, carbon and vitamins sources maltose is a n energy source. Chloramphenicol is a broad spectrum antibiotic that inhibits the majority of gram positive and negative bacteria making the medium selective for *T. vaginalis*.

KIT CONTENT AND COMPOSITION

n. 25 vials Anonymous tubes containing ready to use liquid culture medium
n. 25 Labels

COMPOSITION

Component	Amount/Liter	U.M.
Pancreatic Digest of Casein	12.0	g
Liver extract	2.0	g
Yeast Extract	5.0	g
L-cysteine HCl	1.0	g
Maltose	2.0	g
Agar	1.0	g
Cloramphenicol	0.1	g
Methylen Blue	3.0	mg
Horse Normal Serum	60	ml
pH	6.0 ± 0.2	

Preservation and stability

 = Storage temperature 2-8 °C

Store the reagent at 2-8° C avoiding direct light. Do not freeze. If stored as described above, reagents are stable until the expiration date printed on the label.

Note: During storage some sediments may accumulate on the bottom of the tube (small fibers or flocs) well visible upon shaking. Their presence does not affect the performances of the product

SAMPLE COLLECTION

A) Inoculate a Gynecult tube directly with the patient sample

Women:

- Take the sample with a swab from vagina or the posterior fornix
- Transfer the sample, immediately after the withdrawal of the sample, in the Gynecult tube rotating the swab in the broth
- Remove the swab from the broth

Men:

- Inoculate the Gynecult tube with concentrated urine, urethral swab, a drop of semen or prostatic fluid.
- Shake vigorously

B) Place on the tube the label filled with the personal details of the patient.

C) Incubate the tube at 37°C

D) Maintain the incubation up to 72h

READING

If the sample is positive, after the incubation the broth is cloudy (both with *T. vaginalis* and yeasts). Take a sample from the culture and put it on a glass slide; examine three visual fields with a 40X magnification. At microscope *T. vaginalis* can be seen as a mobile parasite pear-shaped, while yeasts have round or oval cells often showing budding.

T. vaginalis: the presence of a single microorganism, after a 24 hours incubation, is significative of infection.

Yeasts: the result is significative in presence of at least 6 cells. In the presence of a positive result is necessary to proceed with the species identification.

PERFORMANCES CHARACTERISTICS

Culturing is the most sensitive method for *Trichomonas vaginalis* diagnosis. To increase the test sensibility is advisable to perform the direct microscopy identification (wet and/or stained smears can be examined).

Limitation of the procedure

Broth culture does not discriminate between *Trichomonas* and yeasts (further test must be necessary). A positive sample does not discriminate between *Candida* patients with disease and asymptomatic carriers. Using the culturing method is, however, always necessary the the wet smear examination (because the sample may contain non-viable microorganisms).

Quality Control

Positive

C. albicans ATCC 18804
T. vaginalis ATCC 30001

Negative

S. aureus ATCC 25923
E. coli ATCC 25922

Inocula must be observed under the microscopy (40X magnification) before and after incubation, to confirm an increase in the number of the cells. Test of repeatability (within run) and reproducibility (between run) carried out with three different lots of medium with the Quality Control strains, at defined concentrations, have given the expected results on the 98% of cases.

WARNINGS

1. All human samples must be handled and disposed as potentially infectious materials.
2. The use of the not inoculated medium does not present particular risks. After inoculation the medium has to be treated as potentially infectious.
3. The kit should only be used by qualified and properly trained technical personnel.
4. Diagnoses shall only be carried out by authorised and qualified personnel. Interpretation of the test results should be made taking into consideration the patient history, and, if necessary, the results of other tests performed.
5. It is recommended to handle the reagent according to the rules of good laboratory practice and to use appropriate personal protective equipment.
6. Comply with national directives on occupational safety and quality assurance.
7. Use equipment that comply with current standards.
8. Do not use the culture medium after the expiry date

Reporting of serious incidents

Please inform the manufacturer (through your distributor) and the competent authority of the member state of the European Union in which the user and/or patient is established, of cases of serious incident that has occurred in relation to the device. For other jurisdictions, reports of serious incidents must be made in accordance with the regulatory requirements of the home Member State. By reporting serious incidents, you help provide more information about the safety of your in vitro medical diagnostic device.

PRECAUTIONS

Reagent disposal

Disposal must be carried out in accordance with current national, regional or European governmental provisions regarding waste.




Caution in the use of reagents

In addition to risk indications related to the active components, reagents may contain inactive components such as preservatives and detergents. The total concentration of these components is lower than the limits reported in the EC directive 1272/2008 and subsequent amendments and additions. However, It is recommended to handle reagents according to the standards of good laboratory practice

Bibliography

- 1- **Kupferberg AB, Johnson G., Sprince H.** (1948) Nutritional Requirements of *Trichomonas vaginalis*. Proc Soc Exp Biol Med 8, 67:304-308
- 2- **Murray PR, Baron EJ, Pfaller MA, Tenoer FC, Tenover RH.** (1995) Manual of clinical microbiology. 6th ed. American Society for Microbiology; 1995. p. 1213-7.
- 3- **Gelbart SM, Thomason JL, Osypowski PJ, Kellett AV, James JA, Broekhuizen FF.** (1990) Growth of *Trichomonas vaginalis* in commercial culture media. J. Clin Mikrobiol.;28:962-4.
- 4- **Fouts AC, Kraus SJ.** (1980) *Trichomonas vaginalis*: reevaluation of its clinical presentation and laboratory diagnosis. J Infect Dis 1980;141:137-43.
- 5- **Isenberg HD.** (1992) Clinical microbiology procedures hand-book. American Society for Microbiology; 1992. ch. 2:7.9.3.1-6.
- 6- **McMillan A.** (1990) Laboratory diagnostic methods and cryo preservation of *Trichomonads*. In: Honigberg BM, editor. *Trichomonads parasitic in man*. New York: Springer-Verlag;
- 7- **Räsänen S, Eskelinen M, Merilä M, Kaartinen M.** (1982) Diagnosis of *Candida colpitis* using a semiliquid culture medium. Ann Chir Gynaecol;71:340-3

Symbols used for IFU and Packaging

 In vitro diagnostics medical device	 Manufacturer
 Catalog number	 Instruction for Use
 Lot Number	 Storage Temperature
 Expiration Date	 Biological Risk

REVISION	DATE	CHANGES
C	09-2022	Modified for IVDR Compliance Regulation (EU) 2017/746